

Sansin

Hemodialysis Concentrate Solutions



江西三鑫医疗科技股份有限公司

JIANGXI SANXIN MEDTEC CO.,LTD.

地址: 江西南昌小蓝经济技术开发区富山大道999号

Add: No.999 Fushan Road, Xiaolan Economic Development

Zone 330200 Nanchang, Jiangxi P.R. China

电话 (Tel): 86-791-85950275 Fax: 86-791-85988819

邮箱 (E-mail): sales@sanxin-med.com

网址 (Web): www.sanxin-med.com





Hemodialysis Concentrate Solution

- High Purity and easy to dissolve
- Strictly control and high effective dialysis
- Reliable quality, Accurate electroly concentration

Safe quality
Striving for perfection



Technical Data:

Hemodialysis Powder



Type	Model	Capacity	Volume after Dilution
SXG-F	SXG-F-A	1284.5g	5.5L
		1751.6g	7.5L
		2335.2g	10L
		5838.6g	25L
	SXG-F-B	11677.3g	50L
		588g	7.0L
	1008g	12L	
	2520g	30L	
	5040g	60L	

Specifications	Dialysis Powder	SXG-F	SXJ-F	SXS-F
	Dialysis Solution	SXG-Y	SXJ-Y	SXS-Y

Dilution Ratio A: B: RO Water 1:1.225:32.775 1:1.26:32.74 1:1.83:34

Applicable Machine	Fresenius, B. Braun	Nipro, Toray	Grabro, Baxter
--------------------	---------------------	--------------	----------------

Dialysate ion concentration (mmol/L)	Na ⁺	137	138	140
	K	2.0	2.0	2.0
	Mg ²⁺	0.5	0.5	0.5
	Ca ²⁺	1.5	1.5	1.5
	Cl ⁻	109.5	109	107
	AC ⁻	3.0	4.0	4.0
	HC ₃ O ⁻	32.0	31.0	35.0

Powder B for online hemodialysis



Type	Model	Capacity	Volume after Dilution	Package form	Machine
SXG-F	SXG-F-B	700g	8.3L	Cartridge	Gambro Nipro
				Bibag	Fresenius

Customerized formula is available.

